

An Open Source 3D Printable and Arduino-based Alarm System for monitoring Ventilator Airflow

Dr. Alex Maag's Capstone Design Group:

- Karl Frohlich
- Ian McNitt
- Noah Marinaro



Background

Patient demand of medical equipment during the COVID-19 Pandemic has crippled the supply:


- Limited to no access to available ventilators
- Italy, Spain and NYC are experiencing drastic shortages along with rural hospitals
- On March 22, the FDA fast tracked converting sleep apnea machines to ventilators to meet demand
 - CPAP
 - BIPAP



Auburn's RE-INVENT Project


Auburn University successfully re-purposed a standard CPAP machine into a functional emergency ventilator:

- Cost approximately \$700 worth of off-the-shelf components
- Inspiration to expiration ratio: 1:3 - 3:1
- Breaths per minute: 10 - 30 BPM
- Air pressure set: 5 - 25 cmH₂O
- Alarm senses deviations from patient parameters



Problems:


- Aggressive development
- Little to no design refining
- Difficulty validating efficacy
- Unknown component and sensor integrity



Capstone Idea:

Airflow Monitor



- Adds an extra line of safety against malfunctions
- Uses pressure sensor to measure flow rate to validate tidal volume
- Implement an alarm system to verify necessary ventilator-patient lung tidal volume



How Mechanical Ventilation Works


Ventilation for COVID 19 Patients:

- An Acute Respiratory Distress Syndrome (ARDS)
- Mechanical Inhalation: Pre-set Volume or Pressure
- Mechanical Exhalation: Passive Air Removal
 - Focus: Exhale




Ventilator Basics

1. **Tidal Volume**
 - a. The volume of air the lungs will receive
 - b. Set by a medical professional
 - c. What we will be measuring
2. **Inspiratory Phase ('Triggering')**
 - a. Inside the ventilator
3. **Expiratory Phase ('Cycling')**
 - a. Our emphasis
 - b. Outside the ventilator
 - c. Difficult to tamper with



Using a Pitot Tube to Measure Flow Rate



Total Pressure = Static Pressure + Dynamic Pressure

We can use Dynamic Pressure to solve for velocity of the airstream!


$$\text{Dynamic Pressure} = \frac{\rho}{2} V^2$$

Solving for Velocity:

$$V^2 = \frac{2 * (P_T - P_S)}{\rho}$$

$$\begin{aligned} \text{Total Pressure} &= P_T \\ \text{Static Pressure} &= P_S \end{aligned}$$

An Open Source 3D Printable and Arduino-based Alarm System for monitoring Ventilator Airflow



Hardware Overview

MPXV7002DP Differential
Pressure Sensor



Piezoelectric
Speaker



Arduino Nano





Pushbutton








Rocker Switch



4x AA Power
Supply Bracket (6v)



Parts Required


Picture	Part	Price (approx) as of May 2020
	Arduino Nano	\$7
	MPXV7002DP Differential Pressure Sensor	\$16
	Piezoelectric Speaker	\$1
	4x AA Power Supply (6V) + 4 AA Batteries	\$5
	Pushbutton + Rocker Switch + Resistor + Wires + Screws + 3mm diameter rubber pipe	\$11 Total = \$40/unit

Procedure

1. Nurse attaches alarm to ventilator and turns it on.
2. Alarm calibrates for the first ~10 seconds as the patient breathes normally.
3. After 10 seconds, the alarm ensures that the patient continues to breath at a similar flow and interval.
4. If breathing is irregular, the alarm rings until a button is pressed by the nurse to silence it.


Software/Circuit Overview

Modeled and Simulated using Autodesk Tinkercad!



How do you build our Alarm System?

Visit our Public Github Repository!

 Search **Link:** <https://github.com/knf37/COVID19-VentilatorAlarm>

knf37 / COVID19-VentilatorAlarm

Unwatch 1

Star 0

Fork 0

<> Code

Issues 0

Pull requests 0

Actions

Projects 0

Wiki

Security 0

Insights

Settings

A 3D Printable and Arduino-based Alarm System for monitoring Ventilator Airflow

Edit

Manage topics

24 commits

1 branch

0 packages

0 releases

1 contributor

Branch: master

New pull request

Create new file

Upload files

Find file

Clone or download


knf37 Delete Z

Latest commit a052b65 2 minutes ago


3D_Printable_Parts	Create gitkeep	1 hour ago
Arduino_Code	Create gitkeep	1 hour ago
Parts_Required	Rename Required_Parts/gitkeep to Parts_Required/gitkeep	4 minutes ago
Pictures	Add files via upload	10 minutes ago
README.md	Update README.md	7 minutes ago

Future Direction

- Improve alarm system with LED light
 - Green = Normal
 - Yellow = Warning
 - Red = Needs Attention
- Transmitter (ESP8266)
 - Inexpensive (\$7)
 - 24 / 7 monitor capability
 - Enhanced safety



Additional Components



Temperature and Humidity Sensor





Oxygen Sensor

Future Aim

Monitoring the patient's health over time

- i. National Institute of Health ARDS Network protocol:
 - $V_T = (\text{PBW})(6 \text{ mL} / \text{kg PBW})$
- ii. Track the breathing over time
 - Use this data to monitor their health over time



- Average Height By State 2020. (2020, April 6). Retrieved May 8, 2020, from <https://worldpopulationreview.com/states/average-height-by-state/>
- Benson, T. Pitot Tube. Retrieved May 8, 2020, from <https://www.grc.nasa.gov/WWW/k-12/VirtualAero/BottleRocket/airplane/pitot.html>
- Burch, T., & Zabala, M. THE AUBURN RE-INVENT PROJECT. Retrieved May 7, 2020, from <http://eng.auburn.edu/reinvent>
- Khoury-Folkens, C. VENTILATOR SOS. Retrieved May 7, 2020, from <https://www.ventilatorsos.com/>
- Mechanical Ventilation. (2016, July). Retrieved May 8, 2020, from <https://craighospital.org/resources/mechanical-ventilation>
- Miller, S., & Wysocki, M. Tidal volumes for ventilated patients should be determined at the endotracheal tube. Retrieved May 7, 2020, from https://www.hamilton-medical.com/dam/jcr:abe1ac16-8732-4bcd-a209-7c5faa076aa8/flow-sns-position_whitepaper-en-689217.00.pdf
- Silva, P. L., & Rocco, P. (2018). The basics of respiratory mechanics: ventilator-derived parameters. *Annals of translational medicine*, 6(19), 376. <https://doi.org/10.21037/atm.2018.06.06>
- Pressure Measuring - MIT E-Vent: MIT Emergency Ventilator. (2020, April 4). Retrieved May 8, 2020, from <https://e-vent.mit.edu/controls/pressure-measuring/>


References

Questions?



Pitot Tube

Glenn
Research
Center



Bernoulli's Equation :

static pressure + dynamic pressure = total pressure

$$(p_s + r \times \frac{V^2}{2}) = p_t$$

Solve for Velocity :

$$V^2 = \frac{2(p_t - p_s)}{r}$$

Determining the Tidal Volume:

Q = Flow Rate (m^3/s)

V = Volume (m^3)

A = Area (m^2)

d = distance (m)

v = Velocity (m/s)

t = Time (s)

$$Q = V / t = A \cdot d / t = A \cdot v$$

$$\text{Tidal Volume } (V_T) = Q \cdot t$$

Measuring Additional Parameters

- Positive End Expiratory Pressure (PEEP)
 - Low PEEP: 3 - 5 cmH₂O
 - High PEEP: 12+ cmH₂O
- Peak Inspiratory Pressure (PIP)
 - Low PIP < 30 cm H₂O
 - High PIP > 40 cm H₂O
- Driving Pressure: $PIP - PEEP < 15$

Potential Pressure Readings:

- Target PEEP: 15 cmH₂O
- Target PIP: 30 cmH₂O
- With an I/E ratio of 2:1 and target PEEP:
 - Male = 424 mL / Female = 317 mL